







The outer dimensions of the IC should be the same as those of a credit card, ATM card or ID card (85.6 mm × 53.98 mm) with a radius of 2.88–3.48 mm

<p>Front</p>	<div style="display: flex; justify-content: space-between; align-items: center;"> <div data-bbox="581 478 751 579">  <p>Tyber Medical <small>Expand Your Brand</small></p> </div> <div data-bbox="829 478 1109 562">  <p>Implant Card MR Conditional</p> </div> </div> <div style="margin-top: 10px;">  _____  _____  _____ _____ _____ </div>
<p>Back</p>	<div style="display: flex; align-items: flex-start;"> <div data-bbox="548 873 621 940" style="margin-right: 10px;">  </div> <div data-bbox="540 863 1149 1024"> <p>This person is implanted with a Tyber Medical Plating System and can safely undergo an MR exam only under very specific conditions. Scanning under different conditions may result in injury. Full MRI safety information is available in the MRI Safety Information section of the Instructions for Use (IFU), which can be obtained at www.tybermedical.com/eifu.</p> </div> </div> <div style="border: 1px solid black; padding: 10px; margin-top: 10px; text-align: center;"> <p>Place Product Patient Label Here</p> </div> <p style="text-align: right; margin-top: 10px;">LBL-TM202303 Rev A-01</p>